





4 LEDs curing light





Wavelength:380-515nm  
3 Blue,1 Violet,10w LED



Smart controlling chip  
keeps light output consistent  
Light intensity:  
1000mW/cm<sup>2</sup>,1500mW/cm<sup>2</sup>,2300mW/cm<sup>2</sup>  
600 mW/cm<sup>2</sup>(detect mode)



360° rotatable head




Intensity digital display




## Multiple Curing Options

### Ramp mode



### Pulse mode



### Caries detection



# *Curing Pen*

*4 LEDs curing light*



## Technical parameter:

Dimension: 24\*21\*230(mm)

Net Weight: 109g

Light source head size:  $\Phi 13\text{mm}$

Light output lens size:  $\Phi 9.3\text{mm}$

Power supply: >1200 times

2300mW/cm<sup>2</sup>, 3 Sec