




HyperLight

PORTABLE X-RAY UNIT



- Clinical Imagery



•


••

•••


- Sharp, high resolution images with a **0.4** focal spot

Source focal spot

Source focal spot



Subject



Eighteenth

Effect of focal spot size

Engineered for **accurate** images

Designed with
safety in mind





Built-in leaded double
shield to protect the operator
from leakage radiation



Locking system to prevent
unintentional exposure

- Compact size and lightweight can make you operate it by one hand only



- Intuitive user interface

- Smart angulation system

Simplify
your
digital
workflow



Maximize
your
diagnostic
potential





Model : HyperLight

Focal Spot : 0.4

Tube Voltage : 65kV

Tube Current : 2.5mA

Weight : 1.8 kg

Battery : 2500mAh

Exposure Time Range : 0.02-2 s

Line voltage : 100V-240V,
50Hz/60Hz