

EQ-V

Easy and convenient cordless obturator
with innovative cartridge

Your Choice for Warm Gutta Percha Obturation!

The new EQ-V system by META BIOMED
for the reliable, convenient and
precise root canal obturation.
Find out more inside!



USA Version



www.metabiomedshop.com



www.facebook.com/Metabiomed.Endodontics



www.instagram.com/metabiomedus

**META
BIOMED**

The EQ-V – Easy and Quick System

The difference is in the details



Manual System

- Precise extrusion
- Easy maintenance




Replaceable / Rechargeable Battery

- No separate tools needed
- Extra battery allows uninterrupted treatments



Chemical Resistant

- Any disinfectants can be used
- No cracks!



Refillable Cartridges

- Gutta percha fits directly inside the cartridge for clean delivery
- One click lock in



360° Ring Button

- Allows easy access to mandible and maxilla



Safe Heating System

- Auto shut-off ensures patients' safety
- Saves battery life



EQ-V Full Kit




EQ-V Fill




EQ-V Pack

Obturation Sequence -



Clinical Cases:



Before

After

Courtesy of
Dr Jenner O. Argueta

Before


After

Courtesy of
Dr Álvaro Barrios Angulo

Before

After

Courtesy of
Dr Walter Vargas Obando



Contents Type	Fill Unit	Pack Unit	Battery	Single Charger	Dual Charger	Packtip (50/.04, 60/.05) (2ea/1pk)	Fill Needle (25G) (6ea/1pk)	Multi Tool	Cleaning Brush	AC/DC Adapter	Power Cord	GP Bar (100pcs/bottle)	Seal Driver	Plunger Seal (2ea/1pk)
Full Kit	0	0	0 (3ea)		0	0	0	0	0	0	0	0	0	0
Fill Kit	0		0 (2ea)	0			0	0	0	0	0	0	0	0
Pack Kit		0	0 (2ea)	0		0				0	0			

Ordering Information

Item Code	Description	Item Code	Description
EQV033	Full kit	EQV-P38	Pack tip (50/04, 60/05)
EQV034	Fill Kit	EQV-P26	Pack tip (40/03,50/04, 60/05, 60/05L)
EQV032	Pack kit	EQV042	Gutta percha bar
EQV-F35	Refillable cartridge 23G	EQV-F13	Plunger assy
EQV-F36	Refillable cartridge 25G	EQV-F24	Plunger seal (2ea)
EQV-F32	Refillable cartridge 23G, 25G	EQV-F25	Plunger seal driver
EQV-P33	Pack tip (40/03, 50/04)	EQV-C30	Battery
EQV-P34	Pack tip (40/03)	EQV-C20	Single charger
EQV-P35	Pack tip (50/04)	EQV-C10	Twin charger
EQV-P36	Pack tip (60/05)	EQV-S02	Cleaning & Tools Kit (1 Brushes, Multi tool)
EQV-P37	Pack tip (60/05L)	SS009	AC Adapter

Specifications

EQ-V Fill	Heating Time:	35 seconds
	Setting Temperature:	160 or 200°C
EQ-V Pack	Heating Time:	in a second
	Setting Temperature:	180 or 230°C
	Continuous Operation:	4 seconds
Battery	Li-ion:	3.6 V / 1900mAh
	Charging Time:	3 hours
AC/DC Adapter	Input:	100-240V~
	Output:	9.0 V d.c. / 2.0 A
Operating Time	Fill:	Continuous for 1 hour
	Pack:	Continuous for 1.5 hours